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Abbott Laboratories,
North Chicago, IL 60064 USA
Abbott Laboratories, Queenborough,
Kent ME11 5EL United Kingdom
© 2008 Abbott Laboratories. All Rights Reserved.
ART16065 Rev. A 07/08

Table of Contents

Important Information About the

AlphaTRAK™ Blood Glucose Monitoring System 2-3

Important Health-Related Information 3-4

The AlphaTRAK Meter

The AlphaTRAK Blood Glucose Monitoring System Components 5

The AlphaTRAK Meter and How It Works 6-7

The AlphaTRAK Meter Display 8-9

Setting Up the AlphaTRAK Meter 10-12

The AlphaTRAK Test Strip

The AlphaTRAK Test Strip 13

The AlphaTRAK Control Solution

The AlphaTRAK Control Solution 14-15

Performing an AlphaTRAK Control Solution Test 16-17

Using the AlphaTRAK Lancets

Obtain a blood sample from the animal 18

AlphaTRAK Blood Glucose Testing

Performing an AlphaTRAK Blood Glucose Test 19-22

LO and HI Readings 22-23

Meter Memory 24-25

Error Messages 26-28

Battery 29

Replacing the Battery 29-30

Caring for the Meter 30

Troubleshooting 31

System Specifications 32

Symbols 33

References 33

Important Information About The AlphaTRAK™ Blood Glucose Monitoring System

Veterinary medical device for use under the direction of a veterinary professional.

Federal law (USA) restricts this device to use by or on the order of a licensed veterinarian.

How the AlphaTRAK Meter Should Be Used

- UNITS OF MEASURE: In the USA the units will be mg/dL. For outside the USA the units will be mmol/L.
- The AlphaTRAK Meter is intended for use outside the body (in other words, *in vitro* diagnostic use only).
- The AlphaTRAK Meter is used for monitoring blood glucose (also known as *blood sugar*) from samples taken from the dog or cat.
- The AlphaTRAK Meter should only be used with AlphaTRAK Test Strips and AlphaTRAK Control Solution. Using other brands of test strips and control solutions with the AlphaTRAK Meter can produce inaccurate results.

Do use the AlphaTRAK Meter for:

- Testing blood glucose in dogs and cats.
- Testing with fresh whole blood samples.
- Testing with fresh whole capillary samples.
- Testing with fresh whole venous samples.
- Testing with whole blood collected in a syringe or tube with EDTA or heparin anticoagulant if testing within 10 minutes.
- Testing with whole blood collected in a syringe or tube with no anticoagulant, if testing is performed immediately.

Do **not** use the AlphaTRAK Meter for:

- Testing of arterial blood.
- Testing with serum.
- Testing with plasma.
- Diagnosis of diabetes.

WARNING:

The AlphaTRAK Meter contains small parts that may be dangerous if swallowed.

Caution: Please read all of the instructions provided in this *AlphaTRAK Blood Glucose Monitoring System User Guide* and practice the testing procedures before using the AlphaTRAK™ Meter. Blood glucose monitoring should be done with the guidance of a veterinarian/veterinary surgeon.

Important Health-Related Information

- AlphaTRAK Blood Glucose Monitoring System should **NOT** be used to diagnose diabetes.
- Glucose values are affected by stress, meals, health, exercise and medication intake.
- Severe dehydration and excessive water loss may cause false low results. If severe dehydration is suspected, consult a veterinarian/veterinary surgeon immediately.
- The injection of solutions containing dextrose will elevate blood glucose levels.
- Low or high blood glucose readings can indicate a potentially serious medical condition.

General Guidelines are as follows:

Consult a veterinarian/veterinary surgeon for the target glucose values for each animal.

Common symptoms of Hypoglycemia/Hypoglycaemia (low blood sugar):¹

- Dogs: lethargy (abnormal drowsiness), weakness, head tilting, ataxia (shaky movements or unsteady gait), seizures, coma.
- Cats: weakness, ataxia, disorientation, convulsions.

Common symptoms of hyperglycemia/hyperglycaemia (high blood sugar):¹

- Dogs and cats: glycosuria (presence of abnormal amounts of glucose in the urine), polyuria (frequent urination), polydipsia (excessive or abnormal thirst), polyphagia (excessive appetite or eating), weight loss.

If the animal is experiencing symptoms that are not consistent with the blood glucose test result, and all instructions have been followed as described in the *AlphaTRAK™ Blood Glucose Monitoring System User Guide*, contact a veterinarian/veterinary surgeon.

- Hypoglycemia/Hypoglycaemia (low blood sugar) values are:
 - Below 65 mg/dL (3.6 mmol/L) for adult cats and dogs.²
- Hyperglycemia/Hyperglycaemia (high blood sugar) values are:
 - Above 250 mg/dL (13.9 mmol/L) for adult cats and dogs.²

If low or high glucose values are obtained and there are no symptoms of hypoglycemia/hypoglycaemia or hyperglycemia/hyperglycaemia, repeat the test with a new test strip.

If repeated tests indicate low or high glucose levels and all instructions have been followed as described in the *AlphaTRAK Blood Glucose Monitoring System User Guide*, contact a veterinarian/veterinary surgeon.

AlphaTRAK™ Blood Glucose Monitoring System Components

The AlphaTRAK Meter and accessories work together to measure the amount of sugar (glucose) in dog or cat blood. The system includes:

- **AlphaTRAK Meter**
- **AlphaTRAK Blood Glucose Monitoring System User Guide**
 - Contains system information and directions
- **Carrying Case**
 - Use to store and carry meter and other system supplies
- **Registration Card (USA only)**
- **One lithium battery, #2032, 3 volt (replaceable)**

Additional items required to complete testing (sold separately)

- **AlphaTRAK Test Strips (50)**
- **AlphaTRAK Control Solution**
- **AlphaTRAK Lancets (optional)**
- **AlphaTRAK Lancing Device (optional)**

The AlphaTRAK™ Meter and How It Works

Data Port

Data management system is not available for the AlphaTRAK Meter.

Display Screen


Displays the test results and other important information.

m Button

- Powers the meter on and off.
- Moves to different mode settings.
- Changes the test strip code (scrolls backward).

Sample Areas (rectangles)

Apply blood or control solution to one edge of sample area of the strip only.

Sample area of test strip

c Button

- Marks a control solution test.
- Changes the test strip code (scrolls forward), date, time, sound options, and reviews stored readings.

Test Strip Port

Insert the top end of an AlphaTRAK test strip here.

Top

AlphaTRAK Test Strip

Bottom

The AlphaTRAK™ Meter Display

System Check Screen

This screen always appears when the meter is first powered on. If the system check screen does not match the example here, do not use the meter. Please call Customer Care (USA and Canada) or Technical Enquiries (UK).

Memory Mode Symbol

Indicates test results stored in the meter memory.

Battery Symbol

Appears when battery needs replacement.

Control Solution Test Result Symbol

Indicates a test result has been marked as a control solution test.

Test Result Area

Displays test results.

Message Area

Displays date, time, number of tests, and software version.

Check User Guide Symbol

Appears with another symbol to indicate the relevant reference in the *AlphaTRAK Blood Glucose Monitoring System User Guide*.

mg/dL (USA),
mmol/L (outside the USA)

Indicates unit of measurement. Unit is locked according to the country where the meter has been purchased.

Thermometer Symbol

Appears if your meter is out of operating temperature range.

Blood Drop and Test Strip Symbols

These symbols provide the animation indicating when the meter is ready for application of blood or control solution.

Check and Set Test Strip Code Symbol

Indicates when to check and set the test strip code.

14-Day Average Symbol


Shows the average value of glucose readings from the last 14 days.

Sound Level Symbols

Indicates high, low, or off sound volume.

Time Setup Symbol

Appears in Set mode when setting the time format.

Date Setup Symbol

Appears in Set mode when setting the date format.

Setting Up The Meter

• To Enter Setup Mode

Start with the meter off. Press and firmly hold the “m” button until the set display screen appears.

• To Exit at Any Time

Press and hold the “m” button until the meter powers off.

• How to Move Through the Screens

Press the “m” button to move through all the setup screens such as sound level, time, date, and so on.

Press the “c” button to change settings within a setup screen.

Setting the Sound

There are three sound volume options: Off, Low, and High (default setting). Press the “c” button until the desired sound level appears. Press the “m” button to move to “Setting the Time.”

OFF

LOW

HIGH

Setting the Time

HOUR

MINUTE

Setting the Date

DAY

MONTH

YEAR

1. With the hour blinking, press the “c” button until the correct hour appears. If using the 12-hour time format, be sure to press the “c” button until the correct choice between “AM” and “PM” appears.
2. Press the “m” button to move to minute setting.
3. With the minute blinking, press the “c” button until the correct minute appears.
4. Press the “m” button to move to “Setting the Date.”

1. With the day blinking, press the “c” button until the correct day appears.
2. Press the “m” button to move to the month setting.
3. With the month blinking, press the “c” button until the correct month appears.
4. Press the “m” button to move to the year setting.
5. With the year blinking, press the “c” button until the correct year appears.
6. Press the “m” button to move to “Setting the Unit of Measurement.”

Unit of Measurement

1. The units of measurement for the AlphaTRAK™ Meter are milligrams per deciliter (mg/dL) for the USA and millimoles per litre (mmol/L) for outside the USA. The units of measurement are locked and cannot be changed.
2. Press the “m” button to move to “Setting the Time Format.”

Setting the Time Format

- Select the 12-hour or 24-hour time format.
1. Press the “c” button to select the preferred time format. The time will automatically update to the time format selected. If using the 12-hour time format, AM or PM will also appear.
 2. Press the “m” button to move to “Setting the Date Format.”

Setting the Date Format

- Select m-d-y or d-m-y.
1. Press the “c” button to select the preferred date format.
 2. Press the “m” button to move to the “Software Version Display.”


Software Version Display

This is the software version in your AlphaTRAK Meter. It is for display only and cannot be changed.


The picture to the left is an example. The meter may have a different number. To exit the software version display screen, press the “m” button for a few seconds until the meter powers off.

The AlphaTRAK™ Test Strip

Apply blood or control solution to **only one edge of the sample area** of the test strip per test.

Important: Use only one sample area of the test strip per test. Do not apply blood to both sample areas. Test strips may be used only once. Discard used test strips. The test strip can continue to be filled for up to 60 seconds. Be sure to reapply the sample to the same sample area.

WARNING: The test strip cap or vial cap contains drying agents that may be harmful if inhaled or swallowed and may cause skin or eye irritation.

Important Test Strip Information

- The AlphaTRAK Meter has an operating range of 40° – 104° F (5° – 40° C).
- Store the test strip package in a dry place between 39° – 86°F (4° – 30° C). Use test strips only within the system operating temperature range.
- Keep away from direct sunlight and heat.
- Store the test strips in their original vial only; never transfer them to another vial or any other container.
- Never store individual test strips outside the vial.
- After removing an AlphaTRAK Test Strip from the vial, immediately close the vial cap tightly.
- With clean, dry hands, gently touch the test strip when removing it from the vial or inserting it into the meter.
- Do not use test strips beyond the expiration date. This may cause inaccurate results.
- Do not bend, cut, or alter AlphaTRAK Test Strips.

Note: For further information, refer to the AlphaTRAK Test Strip package insert.

The AlphaTRAK™ Control Solution

AlphaTRAK Control Solution is a red liquid that contains a fixed amount of glucose that is used to:

- Practice testing without having to use animal blood.
- Ensure that the meter and test strips are working together properly.

A control solution test should be performed:

- After setting up the AlphaTRAK Meter and before the first blood test.
- When a new vial of test strips is opened.
- When blood sample results are questionable and confirmation that the monitor and test strips are working together properly is required.
- If the test strip vial was not closed properly.
- If the meter is dropped.
- Upon the advice of a veterinary professional.

AlphaTRAK Control Solution is good for three months after opening the bottle or until the expiration date printed on the label, whichever comes first.

Important: Count forward three months from the date a new bottle of control solution is opened. This is the discard date. Write this date on the side of the control solution bottle.

Important Control Solution Information

- Use only AlphaTRAK™ Control Solution with the AlphaTRAK Meter.
- Gently invert the AlphaTRAK Control Solution bottle two to three times before each use.
- To avoid contamination of AlphaTRAK Control Solution, do not touch the tip of container to any surface and discard the first two drops.
- Replace the cap on the control solution bottle immediately after use.
- Do not use control solution past the expiration date.
- Do not add water or any liquid to control solution.
- Control solution tests are specified to be accurate only when used between 59° – 104° F (15° – 40° C).

Caution: Results from AlphaTRAK Control Solution tests do not reflect the animal's blood glucose level.

- The control solution test results should fall within the appropriate range printed on the test strip vial.

Important: The control solution range is a target range for control solution only. It is not a target range for the animal's blood glucose level.

- If control solution results are out of range, repeat the test.
- If the control solution test results continue to fall outside the appropriate range, do not use the AlphaTRAK Meter to test the animal's blood glucose. Contact Customer Care (USA and Canada) call (888) 299-7416 or Technical Enquiries (UK) call + 44 (0) 1628 644354.
- For more information on how to obtain AlphaTRAK Control Solution, contact a veterinarian/veterinary surgeon, Customer Care (USA and Canada) or Technical Enquiries (UK).

Performing an AlphaTRAK Control Solution Test

1

To avoid contaminating the test strip, ensure that hands are clean and dry. Insert the top of the test strip into the meter until it stops. Press the “m” button to power on the meter.

2


The code number on the display must match the code number on the side of the test strip vial.

Note: Be ready to adjust the code on the meter to match the CODE DOG or CODE CAT number on the test strip vial.

Wait for the meter to display the Blood Drop and Test Strip symbols.

3

Sample area of test strip

Gently touch the edge of one test strip sample area to the control solution. The test strip acts like a sponge and pulls the control solution into the strip through the sample area. If the sound prompt is on, the meter will beep when the test strip is full. A short line then moves clockwise around the display screen while the sample is processing.

Important: Do not apply control solution to both sample areas.

16

The result will appear on the display screen, and if the sound is on, the meter will beep twice. To mark the result as a control solution result, press and hold the “c” button for two seconds until the control solution symbol (CONT) appears at the top of the display screen.

If a blood glucose test is marked as a control solution test by mistake, unmark it by pressing and holding the “c” button again for two seconds.

4

on

Compare the control solution test result to the appropriate range printed the side of the test strip vial. The control solution result should fall within the appropriate range.

Out of range test results may be caused by:

- Expired or contaminated control solution.
- Expired or contaminated test strip.
- An error in performing the test.
- Concentrated or diluted control solution.
- Code on test strip vial does not match code set in meter.
- Meter malfunction.
- Control solution test performed outside 59° – 104°F (15° – 40°C).

See page 31 for troubleshooting.

To power off the meter, press and hold the “m” button.

17

Using the AlphaTRAK™ Lancet and obtaining a blood sample from the animal

Consult a veterinary professional for a recommended sampling method and instructions. There are many sites to obtain capillary blood samples: the marginal ear vein in both cats and dogs; paw pad in cats and dogs; leg callus, usually dogs; inner or outer lip, dogs only.

To obtain a capillary sample:

1. Make sure the sample site is clean and dry.
2. Warming the area will help increase blood flow. This can be done by applying a warm cloth to the sample site, or by gently rubbing the area. If using a warm cloth, ensure the cloth is dry as water may dilute the blood sample, or place the cloth inside a plastic bag. Be sure the cloth is not so hot that it burns the animal.
3. Remove the round cover from the lancet and prick the sampling area. (Optional: a thin layer of petroleum jelly may be applied prior to lancing the site. This can help the blood form into a droplet and avoid spilling into the hair. Optional: the AlphaTRAK lancing device may be used to assist in obtaining a blood sample. See instructions for use provided with the AlphaTRAK lancing device.)
4. For test instructions please see page 19.
5. Once the test has been completed, gently apply pressure to the sample area with sterile gauze or a cotton ball/cotton wool to help stop the blood flow.

If you accidentally prick yourself with a lancet, please follow normal first aid procedures. If bleeding persists, please contact your local healthcare professional.

Performing an AlphaTRAK Blood Glucose Test

1

Prepare the AlphaTRAK™ Meter
BEFORE TURNING THE METER ON, insert a new test strip into the meter until it stops. Press the “m” button to power on the meter.

Note: If the test is not started within two minutes, the meter will power off. To restart the meter, press the “m” button to power on the meter.

Setting the Test Strip Code

When the meter powers on, four display screens will appear in the following order:

1. System Check

2. Date and Time

3. Code Number

4. Blood Drop and Test Strip Symbols


Note: To ensure accurate test results, make sure the code number on the meter always matches the number in the outlined box on the test strip vial for either dog or cat.

- When the word “COdE” and a code number is displayed, use the “c” button to scroll forward or the “m” button to scroll backward until the number matches the number on the test strip vial. The AlphaTRAK™ Meter has a total of 50 codes.
- After setting the code number, the Blood Drop and Test Strip symbols will appear on the display screen.
- To return to the code number display screen, press the “c” button.

Caution: Take care to not press the “c” or “m” buttons when holding the meter while performing the blood test. This may change the code and provide an incorrect reading or an Error 6.

2

The test strip acts like a sponge and pulls the blood into the strip through the sample area.

Sample area of test strip

- **Do not** press the test strip against the sample site.
- **Do not** scrape the blood onto the test strip.
- **Do not** apply blood to the test strip when the test strip is not inserted into the meter.
- **Do not** allow blood, control solution or foreign objects to enter into the test strip port of the meter.

Waiting for Result

- Do not pull the test strip away from the sample until *one* beep is heard **or** the short line moving clockwise on the meter screen is observed. This means the sample is adequate and the meter is reading the glucose.
- If after five seconds the meter displays a short dotted line moving back and forth on the display screen, an adequate sample has not entered the test strip. Add more sample to the same area of the test strip within 60 seconds from the time of the first application.

Equivalent to
5.9 mmol/L

- The test result is complete when the blood glucose test result is shown on the display screen. If the sound prompt is on, two beeps will be heard.

Note: To get another test result prior to meter powering off:

- Press “m” to manually power off the meter.
- Remove used test strip.
- Insert new test strip.
- Press the “m” button to power on the meter.
- Follow the steps to perform a blood glucose test.

LO and HI Readings

Important: Low or high blood glucose readings can indicate a potentially serious medical condition. The AlphaTRAK™ Meter displays results from 20 to 750 mg/dL (1.1 to 41.7 mmol/L).

LO Reading

If the test result is lower than 20 mg/dL (1.1 mmol/L), **LO** will appear on the meter display screen.

LO Reading *with Symptoms*

If a **LO** reading is displayed and symptoms such as lethargy (abnormal drowsiness), weakness, head tilting, ataxia (shaky movements or unsteady gait), seizures, or coma in dogs and weakness, ataxia, disorientation, or convulsions in cats¹ are observed, **then** follow a veterinarian/veterinary surgeon's recommendation to treat hypoglycemia/hypoglycaemia (low blood glucose).

LO Reading *without Symptoms*

If a **LO** reading is displayed, but no symptoms of low blood glucose are observed, run a control test to see if meter and test strips are working properly together (reading within range). If

control is within range, **then** retest with a new test strip. If a **LO** reading is again displayed, follow a veterinarian/veterinary surgeon's recommendation to treat hypoglycemia/hypoglycaemia.

HI Reading

If the test is above 750 mg/dL (41.7 mmol/L), **HI** will appear on the display screen. This indicates severe hyperglycemia/hyperglycaemia (high blood glucose).

HI Reading *with Symptoms*

If a **HI** reading is displayed and symptoms such as glycosuria (presence of abnormal amounts of glucose in the urine), polyuria (frequent urination), polydipsia (excessive or abnormal thirst), polyphagia (excessive appetite or eating), or weight loss¹ are observed, **then** follow a veterinarian/veterinary surgeon's recommendation to treat hyperglycemia/hyperglycaemia.

HI Reading *without Symptoms*

If a **HI** reading is displayed, but no symptoms of high blood glucose are observed, run a control to see if meter and test strip are working properly together (reading within range). If control is within range, **then** retest with a new test strip. If a **HI** reading is again displayed, follow a veterinarian/veterinary surgeon's recommendation to treat hyperglycemia/hyperglycaemia.

Meter Memory

Access and Exit Memory

- **To Access**

Start with the meter off. Press the “m” button and release; the meter will display four screens in the following order:

1. System Check Screen
2. Date and Time
3. Code Number
4. Blood Drop and Test Strip Symbol

When the blood drop and test strip symbol is displayed, hold the “m” button and release; the first memory screen to appear will be the 14-day average.

- **To Exit**

Press and hold the “m” button until the meter powers off. Memory mode can be exited at any time

Equivalent to
7.1 mmol/L

14-Day Average

The 14-day average includes all of the readings from the last 14 days starting from the current date and time shown on the meter. (Control solution, LO, and HI readings are not included in the average.)

14d avg = your 14-day average

n = indicates the number of results included in the 14-day average.

Test Readings

The meter stores the last 250 test readings with time and date. Scroll through all of the readings by pressing and holding the “c” button. The first reading is the most recent result.

Equivalent to
4.5 mmol/L

Blood Glucose Test

A result without a control solution symbol (\checkmark CONTR) indicates that it is a blood glucose reading. The date and time of the test will appear on the bottom of the screen.

Equivalent to
5.7 mmol/L

Control Solution Test

A value with (\checkmark CONTR) indicates that this reading is marked as a control solution test. The date and time of the control solution reading will appear on the bottom of the screen.

Equivalent to
5.9 mmol/L

Special Notes

Temperature: A result with a thermometer symbol indicates that the reading was taken outside of the meter's specified operating temperature range and may not be accurate. This value is included in your 14-day average.

LO Reading: A LO reading (under 20 mg/dL / 1.1 mmol/L) is included in the memory as LO.

HI Reading: A HI reading (over 750 mg/dL / 41.7 mmol/L) is included in the memory as HI.

Error Messages

Error 1

Possible Causes

- Sample is too small.
- Problem with the test strip.
- Problem with the meter.
- Very low blood glucose (less than 20 mg/dL / 1.1 mmol/L).

What To Do

1. If symptoms such as lethargy (abnormal drowsiness), weakness, head tilting, ataxia (shaky movements or unsteady gait), seizures, or coma in dogs and weakness, ataxia, disorientation, or convulsions in cats¹ are observed, follow a veterinarian/veterinary surgeon's recommendation for treating hypoglycemia/hypoglycaemia.
2. Conduct a control solution test using a new test strip.
 - A. If the results of the control solution test are **within** the appropriate range printed on the side of the test strip vial, retest using blood and a new test strip.
 - B. If the results of the control solution test are **not within** the appropriate range printed on the side of the test strip vial, or the error persists, call Customer Care (USA and Canada) or Technical Enquiries (UK).

Error 2

Possible Causes

- Problem with the test strip.
- Problem with the meter.
- Very high glucose (above 750 mg/dL / 41.7 mmol/L).

What To Do

1. If symptoms such as glycosuria (presence of abnormal amounts of glucose in the urine), polyuria (frequent urination), polydipsia (excessive or abnormal thirst), polyphagia (excessive appetite or eating), or weight loss¹ are observed, follow a veterinarian/veterinary surgeon's recommendation for treating hyperglycemia/hyperglycaemia.
2. Conduct a control solution test using a new test strip.
 - A. If the results of the control solution test are **within** the appropriate range printed on the side of the test strip vial, retest using blood and a new test strip.
 - B. If the results of the control solution test are **not within** the appropriate range printed on the side of the test strip vial, or the error persists, call Customer Care (USA and Canada) or Technical Enquiries (UK).

Error 3

Possible Causes

- Incorrect test procedure. For example, putting blood on the test strip before inserting the test strip into the meter or applying blood before the Blood Drop and Test Strip symbols display.
- Problem with the test strip.
- Problem with the meter.

What To Do

1. Ensure that the Blood Drop and Test Strip symbols are displayed on the screen before applying blood or control solution.
2. Conduct a control solution test using a new test strip.
 - A. If the results of the control solution test are **within** the appropriate range printed on the side of the test strip vial, retest using blood and a new test strip.
 - B. If the retest results of the control solution test are **not within** the appropriate range printed on the side of the test strip vial, or the error persists, call Customer Care (USA and Canada) or Technical Enquiries (UK).

Error 4

Possible Causes

- Problem with the test strip.
- Problem with the meter.

What To Do

1. Conduct a control solution test using a new test strip.
 - A. If the result of the control solution test is **within** the appropriate range printed on the side of your test strip vial, retest using blood and a new test strip.
 - B. If the results of the control solution test are **not within** the appropriate range printed on the side of the test strip vial, or the error persists, call Customer Care (USA and Canada) or Technical Enquiries (UK).

Note: Error 5 is not used.

Error 6

Possible Cause

A button on the meter was pushed while applying blood sample or control solution.

What To Do

1. Retest using a new test strip. Do not touch the “m” or the “c” button while applying blood or control solution to the test strip.
2. If the error persists, call Customer Care (USA and Canada) or Technical Enquiries (UK).

Battery

The AlphaTRAK™ Meter comes with one pre-installed #2032, 3 volt lithium battery. The battery provides enough power to perform about 1,000 tests. When the battery runs low, the battery symbol () appears on every display screen until the battery is changed.

Important: When this symbol appears, replace the battery immediately.

- To avoid losing the time and date settings, replace the battery within one minute of removing the old battery.
- If it takes more than one minute to replace the battery, time and date may require resetting to get the correct 14-day average. If the time and date settings are lost and not reset, new test results will be stored in memory with a time of 12:00 a.m. (00:00) and a date of 00-00.
- Removing the battery does not affect the meter's memory log or previously-stored user settings (such as calibration code, and date format).
- If the time and date settings are lost, a blood glucose test can be performed and results will not be affected.


Replacing the Battery

1

Turn the meter over. Open latch and pull battery door to open. Lift the tab and the battery will come out automatically.

2

Install the new battery with the positive (+) sign visible, facing outward from the meter. Ensure battery is installed on top of tab.

3

To close, insert the battery door tab in place, then press down firmly on latch until it snaps securely.

Note: Batteries need to be properly disposed of. For this purpose, all shops selling batteries provide battery collection boxes on their premises.

Caring for The Meter

Avoid getting dirt, dust, blood, control solution, water, or any other substance in the meter's test strip and data ports.

Clean the outside of the meter using a cloth dampened with one of the following:

- Mild detergent/mild soap and water,
- 70% Isopropyl alcohol (available in the USA), or
- A mix of one part household bleach, nine parts water.

Important: *Never immerse the meter in water or any other liquid solution.*

30

Troubleshooting

Meter does not power on after inserting a test strip.

Probable Cause	What to Do
The battery is dead.	Replace the battery and reset the date and time, if necessary.
The battery is installed incorrectly or there is no battery in the meter.	Check that the battery is installed correctly with the positive (+) sign facing outward from the meter.
Defective meter or test strips.	Call Customer Care in the USA and Canada at (888) 299-7416. For Technical Enquires in the UK call + 44 (0) 1628 644354.
Blood or foreign objects put into test strip port.	Call Customer Care in the USA and Canada at (888) 299-7416. For Technical Enquires in the UK call + 44 (0) 1628 644354.

Test does not start after applying the blood sample.

Probable Cause	What to Do
Blood sample too small.	Add blood to the same sample area of the test strip for up to 60 seconds from the time of the first application. If necessary, repeat the test using a new test strip and a larger blood sample.
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic shutoff (two minutes).	Repeat the test using a new test strip. Wait until the Blood Drop and Test Strip symbols appear on the display screen before applying the blood sample.
Defective meter or test strips.	Call Customer Care in the USA and Canada at (888) 299-7416. For Technical Enquires in the UK call + 44 (0) 1628 644354.

31

System Specifications

Assay Method	Coulometric electrochemical sensor
Automatic Shutoff	Two minutes after last user action
Battery Life	~1,000 tests
Calibration	Plasma equivalent
Hematocrit/Haematocrit	15% to 65%
Measurement Units	mg/dL for the USA and mmol/L for outside the USA
Memory	250 blood glucose and control solution tests with date and time
Operating Relative Humidity	5% to 90%
Operating Temperature	40° – 104° F (5° – 40° C)
Power Source	One lithium battery, #2032, 3 volt (replaceable)
Result Range	20 to 750 mg/dL (1.1 to 41.7 mmol/L)
Sample	Whole blood, capillary or venous
Sample Size	0.3 microliter (300 nanoliters)
Size	9.7 cm x 5.3 cm x 2.5 cm 3.8 in x 2.1 in x 1.0 in
Test Time	Average of 15 seconds
Weight	2.2 oz. (62.6 g) including battery

Symbols

	Consult instructions for use		Manufactured by
	Temperature limitation	LOT	Batch code
	Do not re-use		Use by Expressed as CCYY-MM or CCYY-MM-DD
IVD	<i>in vitro</i> diagnostic medical device	SN	Serial number
STERILE R	Sterilized using radiation	REF	Catalog number
	Caution refer to accompanying document		

- ¹ BSAVA *Manual of Canine and Feline Endocrinology*, Third Edition, British Small Animal Veterinary Association. Edited by Carmel Mooney and Mark E Peterson. 2004.
- ² Ettinger SJ, Feldman EC. *Textbook of Veterinary Internal Medicine*. Vol 2. WB Saunders Company. pp 1529-1534. 1995

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